

## cECL Western Blot Kit

**Product Number: CWB01**

---

### Shipping and Storage

Store at 2-8 °C in dark.

### Components

Component	CWB01	CWB01
cECL-A(Luminol)	25 ml	125 ml
cECL-B(Peroxide)	25 ml	125 ml

### Description

The cECL Western Blot Kit is a low background substrate detection kit used in immunoblotting experiments in conjunction with horseradish peroxidase (HRP). This product can undergo chemical reactions and emit light under the catalysis of HRP, and can be used to detect biomolecules such as proteins fixed on membranes. Its high sensitivity can detect ng level antigens, and the luminescent signal is strong and persistent. It can be detected using instruments such as X-ray film exposure or chemiluminescence imaging.

### Note

1. During contact with the membrane, please wear gloves and use clean equipment such as tweezers to avoid protein contamination and high background.
2. Under dark conditions, the prepared chemiluminescence detection substrate working solution can be stably stored at room temperature for 8 hours. Sunshine or other strong light can affect the working fluid, so prolonged exposure to strong light should be avoided. Short term exposure to normal laboratory lighting does not affect the use of working fluids.

### Protocol

1. After the second antibody incubation is completed, wash the imprinting film thoroughly.
2. According to the required amount, mix cECL-A and cECL-B in a 1:1 ratio and equal volume to prepare a luminescent detection substrate working solution (approximately 1 ml of working solution is used for an 8 cm x 6 cm membrane).
3. Discard the washing buffer and drop the luminescent substrate working solution onto the imprinting film, ensuring that the working solution covers the entire film. Incubate at room temperature for 3-5 minutes.
4. Use a pipette to remove excess luminescent substrate working solution and place the imprint film between two clean plastic films. This process should be completed carefully to avoid the formation of bubbles between films.
5. Expose X-ray film in a darkroom or place the film in a chemiluminescence imager and perform testing according to the instrument manual.

### Schedule

Problem	Reason	Resolvent
Film inversion (white stripes, black background)	Excessive HRP in the system	Dilute HRP markers at least 10 times or more
Brown or yellow stripes appear on the membrane		
Strong luminescence seen in the darkroom		
The duration of the luminous signal is		

too short		
Weak or no signal	Excessive HRP in the luminescent reaction system leads to rapid substrate consumption, resulting in rapid signal reduction	Dilute HRP markers at least 10 times or more
	Insufficient antigen/antibody levels	Increase antigen/antibody usage
	Low protein transfer rate	Optimize transfer system
High background	Excessive HRP in the system	Dilute HRP markers at least 10 times or more
	Insufficient closure	Optimize closed programs
	Improper selection of sealing reagents	Choose another blocking reagent
	Insufficient flushing	Increase flushing time and frequency
	Overexposure	Reduce exposure time
	Antigen/antibody concentration too high	Reduce antigen/antibody usage concentration
The protein bands are punctate	Protein transfer failure	Optimize the transfer process
	Membrane imbalance	Handle the film according to the instructions
	There are bubbles between the film and film	Remove all bubbles before exposure
Non specific bands appear (high background, short signal maintenance time)	There are too many HRPs in the system	Dilute HRP markers
Non specific bands appear (with a clean background and normal signal maintenance time)	Excessive dosage of primary antibody	Further dilution of primary antibody
	SDS leads to non-specific binding	Avoid using SDS during the experiment