

# Avian Leukosis Virus (ALV) Universal Nucleic Acid Detection Kit (Fluorescent PCR Method)

Product Number:DTK021

---

## Shipping and Storage

1.  $-20^{\circ}\text{C}\pm 5^{\circ}\text{C}$ , stored in the dark, transported, and subjected to repeated freeze-thaw cycles no more than 5 times, with a validity period of 12 months.
2. The collected or processed samples should be stored at  $2^{\circ}\text{C}\sim 8^{\circ}\text{C}$  for no more than 24 hours; If long-term storage is required, it should be stored at  $-70^{\circ}\text{C}$  or below, with no more than 3 freeze-thaw cycles.

## Component

Component	50T
ALV reaction solution	500 $\mu\text{L}\times 2$
Enzyme solution	50 $\mu\text{L}$
ALV positive quality control product	50 $\mu\text{L}$
Negative quality control product	250 $\mu\text{L}$

## Description

This kit uses a pair of avian leukemia virus specific primers, combined with a specific fluorescent probe, to perform in vitro amplification and detection of avian leukemia virus RNA using one-step fluorescent RT-PCR technology, for clinical pathogen diagnosis of suspected infected individuals.

## Application

Avian Leukosis Virus (ALV) is an important pathogen that causes tumor diseases in poultry. It can also lead to severe cellular and humoral immune suppression as well as growth inhibition in infected organisms, making it prone to co infection or co infection with other pathogens such as Marek's disease virus (MDV) and infectious bursal disease virus (IBDV). At present, there are no effective vaccines and drugs for the prevention and treatment of ALV, mainly through continuous testing and elimination of positive chickens to continuously purify the chicken population. Fluorescent PCR technology has become an important means of detecting avian leukemia pathogens due to its high sensitivity, strong specificity, and easy operation.

This kit is suitable for detecting avian leukemia virus in whole blood, plasma or serum with white blood cells, spleen, liver, kidney and other specimens, and is suitable for auxiliary diagnosis of avian leukemia virus infection.

## Applicable instruments

ABI, Agilent MX3000P/3005P, LightCycler, Bio Rad, Eppendorf and other series of fluorescence quantitative PCR detectors.

## Sample collection

Suspected diseased chicken whole blood, plasma or serum with white blood cells; The spleen, liver, and kidneys of suspected sick chickens.

## Protocol

### 1. Sample processing (sample processing area)

#### 1.1. Sample Preparation

Weigh approximately 1g of each tissue from three different locations, cut and mix it with surgical scissors, and then take 0.5g and grind it in a grinder. Add 1.5mL of physiological saline and continue grinding. After homogenization, transfer it

**For Research Use Only**

to a 1.5mL sterile centrifuge tube and centrifuge at 8000rpm for 2 minutes. Take 100µL of the supernatant and transfer it to a 1.5mL sterile centrifuge tube; Take 100µL of serum from the whole blood sample after coagulation and place it in a 1.5mL centrifuge tube.

**1.2. Nucleic acid extraction**

We recommend using our company's nucleic acid extraction or purification reagents (magnetic bead method or centrifugal column method) for nucleic acid extraction. Please follow the reagent instructions for operation.

**2. Reagent preparation (reagent preparation area)**

Based on the total number of samples to be tested, the required number of PCR reaction tubes is N (N=number of samples+1 negative control tube+1 positive control tube); For every 10 parts of the reaction tube, an additional 1 part is prepared. The preparation of each test reaction system is shown in the table below:

reagent	ALV reaction solution	Enzyme solution
Dosage (sample size N)	19µL	1µL

Transfer the mixed test reaction solution into a PCR reaction tube at a concentration of 20uL per tube.

**3. Sample addition (sample processing area)**

Take 5µL of the nucleic acid, positive control sample, and negative control sample extracted in step 1, and add them to the corresponding reaction tubes. Cover the tubes, mix well, and briefly centrifuge.

**4. PCR amplification (nucleic acid amplification zone)**

4.1. Place the reaction tube to be tested in the reaction tank of the fluorescence quantitative PCR instrument;

4.2. Set the channel and sample information, and set the reaction system to 25µL;

Fluorescence channel selection: Detection channel (Reporter Dye) FAM, Quencher Dye NONE, please do not select ROX reference fluorescence for ABI series instruments, select None.

4.3. Recommended loop parameter settings:

step	Cycles	Temperature	Time	Collect fluorescence signals
1	1 cycle	50°C	10min	No
2	1 cycle	95°C	2min	No
3	45 cycles	95°C	15sec	No
		60°C	30sec	Yes

**5. Result analysis and judgment**

**5.1. Result Analysis Condition Setting**

(Please refer to the user manuals of each instrument for setting up, taking the ABI7500 instrument as an example)

After the reaction is complete, the results will be automatically saved. Based on the analyzed image, adjust the Start value, End value, and Threshold value of the baseline (users can adjust them according to their actual situation, with Start value set between 3-15 and End value set between 5-20, so that the threshold line is in the exponential period of the amplification curve, and the amplification curve of negative quality control products is flat or below the threshold line). Click Analyze to automatically obtain the analysis results.

**5.2. Result judgment**

Positive: The Ct value of the detection channel is ≤ 40, and the curve shows a significant exponential growth curve;

Negative: The Ct value of the sample test result is greater than 40 or there is no Ct value.

**Quality control standards**

Negative quality control product: no specific amplification curve or Ct value display;

Positive quality control product: The amplification curve shows a significant exponential growth period, and the Ct value is ≤32;

The above conditions should be met simultaneously, otherwise the experiment will be considered invalid.

**Limitations of detection methods**

1. The results of sample testing are related to the quality of sample collection, processing, transportation, and preservation;
2. Failure to control cross contamination during sample extraction can result in false positive results;
3. Leakage of positive controls and amplification products can lead to false positive results;



## **MEBEP TECH(HK) Co., Limited**

Email: [sales@mebep.com](mailto:sales@mebep.com) Website: [www.mebep.com](http://www.mebep.com)

Tel: +86-755-86134126 WhatsApp/Facebook/Twitter: +86-189-22896756

4. During the epidemic, genetic mutations and recombination of pathogens can lead to false negative results;
5. Different extraction methods have differences in extraction efficiency, which can lead to false negative results;
6. Improper transportation, storage, or inaccurate preparation of reagents can lead to a decrease in reagent detection efficiency, resulting in false negatives or inaccurate quantitative testing results;
7. The test results are for reference only. If a diagnosis is required, please combine clinical symptoms and other testing methods.

### **Note**

1. All operations must be strictly carried out in accordance with the instructions;
2. The various components in the reagent kit should be naturally melted, completely mixed, and briefly centrifuged before use;
3. The reaction solution should be stored away from light;
4. Try to avoid the presence of bubbles during the reaction, and cover the tube tightly;
5. Use disposable suction tips, disposable gloves, and specialized work clothes for each area;
6. Sample processing, reagent preparation, and sample addition should be carried out in different areas to avoid cross contamination;
7. After the experiment is completed, treat the workbench and pipette with 10% hypochlorous acid, 75% alcohol, or a UV lamp;
8. All items in the reagent kit should be treated as contaminants and handled in accordance with the "Biosafety Guidelines for Microbial Biomedical Laboratories".