

Clonorchis Sinensis Detection Kit (Real-Time PCR Method)

Product Number:DTK303

Shipping and Storage

1. -20°C± 5°C, stored in the dark, transported, and subjected to repeated freeze-thaw cycles no more than 5 times, with a validity period of 12 months.
2. The collected or processed samples should be stored at 2 °C~8 °C for no more than 24 hours; If long-term storage is required, it should be stored at -70 °C or below, with no more than 3 freeze-thaw cycles.

Component

Component	50T
CS reaction solution	500µL×2
Enzyme solution	50µL
CS positive quality control product	50µL
Negative quality control product	250µL

Note: Different batches of reagents cannot be mixed.

Description

This kit uses a specific primer for the nucleic acid of Clonorchis sinensis, combined with a specific probe, to amplify and detect the DNA of Clonorchis sinensis in vitro using fluorescence PCR technology. It is used for pathogen diagnosis of suspected infectious materials in clinical practice.

Application

Clonorchiasis is a foodborne zoonotic parasitic disease caused by infection with Clonorchis sinensis (Cs), which seriously endangers human health. The larvae of Clonorchis sinensis need to pass through two intermediate hosts. The first intermediate host includes various freshwater snails, and the second intermediate host includes various freshwater fish and shrimp.

This kit is suitable for detecting Clonorchis sinensis in stool and other samples, and is used for auxiliary diagnosis of Clonorchis sinensis infection. The test results are for reference only.

Applicable instruments

ABI7500, Agilent MX3000P/3005P, LightCycler, Bio Rad, Eppendorf and other series of fluorescence quantitative PCR detectors.

Specimen collection

For fish with clinical symptoms, such as those with a body length of no more than 4cm, take the whole fish. For fish with a body length of 4-6cm, take the internal organs (including kidneys); Fish with a body length greater than 6cm should have their brain, liver, kidney, and spleen taken; For asymptomatic fish, kidney, spleen, cheek, and brain tissues should be taken, and mature female fish also need to take ovarian fluid. Other mammals take 0.5-1mL of water samples.

Protocol

1. Sample processing (sample processing area)

1.1. Sample Preparation

Weigh approximately 1g of each tissue from three different locations, cut and mix it with surgical scissors, and then take 0.5g and grind it in a grinder. Add 1.5mL of physiological saline and continue grinding. After homogenization, transfer it to a 1.5mL sterile centrifuge tube and centrifuge at 8000rpm for 2 minutes. Take 100 µ L of the supernatant and transfer it

For Research Use Only

to a 1.5mL sterile centrifuge tube; Centrifuge the water sample at 8000rpm for 2 minutes, take 100 μ L of the supernatant and transfer it to a 1.5mL sterilized centrifuge tube.

1.2. Nucleic acid extraction

We recommend using our nucleic acid extraction or purification reagents (magnetic bead method or centrifugal column method) for nucleic acid extraction. Please follow the instructions in the reagent manual.

2. Reagent preparation (reagent preparation area)

Based on the total number of samples to be tested, the required number of PCR reaction tubes is N (N=number of samples+1 negative control tube+1 positive control tube); For every 7 samples, an additional 1 sample is prepared. The preparation of each test reaction system is shown in the following table:

reagent	CS Reaction solution	Enzyme solution
Dosage (sample size N)	19μL	1μL

Transfer the mixed test reaction solution into a PCR reaction tube at a concentration of 20uL per tube.

3. Sample addition (sample processing area)

Take 5μL of the nucleic acid, positive control sample, and negative control sample extracted in step 1, and add them to the corresponding reaction tubes. Cover the tubes, mix well, and briefly centrifuge.

4. PCR amplification (nucleic acid amplification zone)

4.1. Place the reaction tube to be tested in the reaction tank of the fluorescence quantitative PCR instrument;

4.2. Set the channel and sample information, and set the reaction system to 25μL;

Fluorescence channel selection: Detection channel (Reporter Dye) FAM, Quencher Dye NONE, please do not select ROX reference fluorescence for ABI series instruments, select None.

4.3. Recommended loop parameter settings:

step	Cycles	Temperature	Time	Collect fluorescence signals
1	1 cycle	95°C	2min	No
2	45 cycles	95°C	15sec	No
		60°C	30sec	Yes

5. Result analysis and judgment

5.1. Result Analysis Condition Setting

(Please refer to the user manuals of each instrument for setting up, taking the ABI7500 instrument as an example)

After the reaction is complete, the results will be automatically saved. Based on the analyzed image, adjust the Start value, End value, and Threshold value of the baseline (users can adjust them according to their actual situation, with Start value set between 3-15 and End value set between 5-20, so that the threshold line is in the exponential period of the amplification curve, and the amplification curve of negative quality control products is flat or below the threshold line). Click Analyze to automatically obtain the analysis results.

5.2. Result judgment

Positive: The Ct value of the detection channel is ≤ 40, and the curve shows a significant exponential growth curve;

Negative: The Ct value of the sample test result is >40 or there is no Ct value.

Quality control standards

Negative quality control product: no specific amplification curve or Ct value display;

Positive quality control product: The amplification curve shows a significant exponential growth period, and the Ct value is ≤32;

The above conditions should be met simultaneously, otherwise the experiment will be considered invalid.

Limitations of detection methods

1. The results of sample testing are related to the quality of sample collection, processing, transportation, and preservation;
2. Failure to control cross contamination during sample extraction can result in false positive results;
3. Leakage of positive controls and amplification products can lead to false positive results;
4. During the epidemic, genetic mutations and recombination of pathogens can lead to false negative results;
5. Different extraction methods have differences in extraction efficiency, which can lead to false negative results;



MEBEP TECH(HK) Co., Limited

Email: sales@mebep.com Website: www.mebep.com

Tel: +86-755-86134126 WhatsApp/Facebook/Twitter: +86-189-22896756

6. Improper transportation, storage, or inaccurate preparation of reagents can lead to a decrease in reagent detection efficiency, resulting in false negatives or inaccurate quantitative testing results;
7. The test results are for reference only. If a diagnosis is required, please combine clinical symptoms and other testing methods.

Note

1. All operations must be strictly carried out in accordance with the instructions;
2. The various components in the reagent kit should be naturally melted, completely mixed, and briefly centrifuged before use;
3. The reaction solution should be stored away from light;
4. Try to avoid the presence of bubbles during the reaction, and cover the tube tightly;
5. Use disposable suction tips, disposable gloves, and specialized work clothes for each area;
6. Sample processing, reagent preparation, and sample addition should be carried out in different areas to avoid cross contamination;
7. After the experiment is completed, treat the workbench and pipette with 10% hypochlorous acid, 75% alcohol, or a UV lamp;
8. All items in the reagent kit should be treated as contaminants and handled in accordance with the "Biosafety Guidelines for Microbial Biomedical Laboratories"