



## Fast Pfu DNA Polymerase

Product Number:PC0302

---

### Storage Conditions

-20°C

### Shelf Life

3 years

### Concentration

5U/μl

### Description

The enzyme is a heat-resistant DNA polymerase with high fidelity, and its fidelity is 8~10 times that of Taq DNA polymerase. The enzyme has not only 5'-3'DNA polymerase activity, but also 3'-5'exonuclease proofreading activity. It can timely remove the base mismatched nucleotides in the DNA synthesis chain, greatly increasing the correct rate of base pairing and ensuring the high fidelity of DNA synthesis. The enzyme has no 5'-3'exonuclease activity. The synthesis speed is 15 seconds 1KB, and the effective amplification length is more than 8KB.

### Source

E. Coli recombinant strain, containing the gene cloned from Pyrococcus furiosus and modified.

### Unit Definition of activity

Take calf thymus DNA as template, react at 72°C for 30min, and the amount of enzyme that can catalyze the polymerization of 10nmol dNTP into polynucleotide fragments is defined as 1 active unit (U)

### Storage buffer

[PH8.0] 20 mM Tris-HCl, 100 mM KCl, 0.1 mM EDTA, 0.1 % NP-40 (v/v), 0.1 % Tween-20(v/v), 1 mM DTT, 50 % glycerol(v/v).

### Reaction buffer

[10×Pfu Buffer (Mg<sup>2+</sup>)] 200 mM Tris-HCl (PH 8.8), 100 mM KCl, 160 mM (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>, 1 % Triton X-100, 1 mg/ml BSA, 20 mM MgSO<sub>4</sub>

### Protocol

User prepared reagent: dNTP mixture; Water (without nuclease); Primers (upstream and downstream); Template

#### 1. Establishment of PCR reaction system

After all the solutions are completely melted, shake gently and fully;

Add into thin-walled PCR tube and operate on ice; fifty μ L the reaction system is as follows:

| Component                        | Addition volume (μl) | Final concentration |
|----------------------------------|----------------------|---------------------|
| 10× Buffer with Mg <sup>2+</sup> | 5                    | 1×                  |
| dNTP Mix,10mM                    | 1                    | 0.2 mM              |
| Upstream primer                  | X                    | 0.1-1.0μM           |
| Downstream primer                | Y                    | 0.1-1.0μM           |
| pfu (5U/μl)                      | 0.5                  | 2.5U                |

**For Research Use Only**



# MEBEP TECH(HK) Co., Limited

Email: [sales@mebep.com](mailto:sales@mebep.com) Website: [www.mebep.com](http://www.mebep.com)

Tel: +86-755-86134126 WhatsApp/Facebook/Twitter: +86-189-22896756

|              |      |             |
|--------------|------|-------------|
| DNA Template | Z    | <0.5µg/50µl |
| Add water to | 50µl | -           |

**Note:** PFU enzyme should be added after dNTPs is added, otherwise the fidelity activity of PFU enzyme may degrade the primer, resulting in the production of non-specific bands and the decline of product yield.

## 2. PCR reaction cycle setting

| Program          | Temperature | Time       | cycles  |
|------------------|-------------|------------|---------|
| Pre denaturation | 95°C        | 2 min      | 1       |
| Denaturation     | 95°C        | 2 min      | } 25-35 |
| Annealing        | 42-65°C     | 0.5-1 min  |         |
| Extend           | 72°C        | 10-15s /kb |         |
| Final extension  | 72°C        | 5 min      | 1       |

## 3. Result test

Take 5 after reaction µ L reaction products were detected by agarose gel electrophoresis.