

## ThermoPlus M-MLV Reverse Transcriptase (RNase H-)

Product Number: RT5060

---

### Shipping and Storage

Store at -20°C.

### Components

Component	RT5060 200KU
ThermoPlus M-MLV Reverse Transcriptase (RNase H-)	1 ml
5×RT Buffer	1 ml

### Description

ThermoPlus M-MLV Reverse Transcriptase(RNase H-), encoded by Moloney Murine Leukemia Virus (M-MLV RT) is an RNA-dependent DNA polymerase that synthesizes the complementary cDNA first strand from a single-stranded RNA template to which a primer has been hybridized. MMLV RT will also extend primers hybridized to single-stranded DNA. Second strand cDNA synthesis can be achieved from some mRNA templates without an additional DNA polymerase. ThermoPlus MMLV is a mixture of several temperature-resistant reverse transcriptase enzymes with good temperature adaptability. The optimal temperature range is 50-60 °C, adapting to a variety of different reaction conditions. It has excellent reverse transcription performance.

### Source

Recombination of E.coli containing Moloney murine leukemia virus reverse transcriptase gene from clone of Moloney murine.

### Concentration

200U/μl

### Features

1. Lack RNase H activity: Weak RNaseH activity High cDNA yield, can get more full length cDNA.
2. Thermal stable: the optimum reaction temperature is 50°C, the highest is 60°C. Can overcome the template RNA secondary structure, and finish the reverse transcriptase experiment smoothly.
3. Wide temperature range: can reverse transcript from 37-60°C, with more than 80% of the highest activity at 42°C-55°C customer can choose the reaction temperature freely.
4. Strong amplification activity: Gene mutation enhanced the binding capacity of the enzyme and RNA. Increased the amplification speed, can obtain the quality cDNA, suitable for cDNA library construction.

### Application

The first-strand cDNA synthesis; RT-PCR.

### Unit definition

One unit of MMLV RT catalyzes the incorporation of 1 nmol of dTTP into acid-insoluble material in 10 minutes at 37°C using oligo(dT)<sub>12-18</sub>-primed poly(A)<sub>n</sub> as a template.

### Storage buffer

20 mM Tris-HCl (pH 7.5), 200 mM NaCl, 0.25 mM EDTA, 0.01% NP-40(v/v), 2.5 mM DTT, 50% glycerol (v/v).

5×RT Buffer: 250mM Tris-HCl (pH 8.3), 15mM MgCl<sub>2</sub>, 375 mM KCl, 50mM DTT.

**For Research Use Only**

**Recommended Reaction Conditions:****The first-strand cDNA synthesis**

1. Add the following reagents to a RNase free PCR tube at room temperature add the MMLV RT last.

Component	Volume
Oligo dT12-18 (1µg/µl) or random primer (50-250ng)	1µl
Total RNA (10ng-5µg) or mRNA(1-500ng)	xµl
dNTP (10mM each)	1µl
DEPC ddH <sub>2</sub> O	(14-x)µl

2. Gently mix and incubate 10 Min at 70°C then chill on ice for 2-10min.
3. Centrifuge for a few seconds then Put the tube into ice and add the next composition :

Component	Volume
5×RT Buffer	4µl
RNasin (40U/µl)	1µl

4. Gently mix and incubate at 50°C for 2 Min (Oligo dT12-18 or sequence especially primer) or at 25°C for 10 min for the random primer.
5. Centrifuge for a few seconds. Add 1µl M-MLV RT (200U/µl) Incubate at 50°C for 50min.
6. Inactivate at 70°C for 10min then get the cDNA.