

## PCRD lateral flow strips

Product Number: PCRD01

---

### Packing

10T/bag

### Description

Used for detecting two genes (or two loci). This test strip is suitable for use with biotin FITC/FAM. The detection of nucleic acid amplification products with third markers is suitable for PCR amplification, isothermal amplification (RPA, LAMP, etc.).

### Protocol

This product uses colloidal gold lateral chromatography double antibody sandwich method to rapidly detect nucleic acid amplification products. The detected nucleic acid amplification product must be labeled with biotin, FITC/FAM, or a third marker at the same time. When the double labeled amplification product flows through the test strip, it can be recognized by colloidal gold labeled antibodies and nitrocellulose membrane immobilized antibodies, forming a red band in the corresponding area.

Structure diagram

